




A wide variety of cross connectors

- Cross connectors are adjustable with pivoting heads
- Ergonomically designed multifunctional instruments
- Connects seamlessly to the thoracolumbar region with the Coral™ Spinal System
- Low profile screw and hook designs
- Occipital plate allows for adjustable rod placement

8mm-28mm (2mm increments)
30mm - 50mm (5mm increments)*

8mm-28mm (2mm increments)

18mm-38mm (2mm increments)

Diameters	Dual Diameter	Lengths
3.5mm	Transition Rod 3.5mm/5.5mm	80mm, 120mm, 240mm (straight) 80mm (Pre-Curved)* 170mm, 205mm*, 240mm* (Pre-Bent Occipital Rod) 385mm (163mm/10mm/212mm), 600mm (203mm/10mm/387mm)*

